


Product Datasheet

Human Kallikrein 5/KLK5 ELISA Kit (orb865336)

Description	Human Kallikrein 5/KLK5 ELISA Kit
Reactivity	Human
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: V23-S293
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA), human milk and saliva
Uniprot ID	Q9Y337
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
KLK5
PicoKine
ELISA Kit
standard
C...