


Product Datasheet

Human ACHE ELISA Kit (orb865337)

Description	Human ACHE ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: E32-L614
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Uniprot ID	P22303
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is around 500 pg/ml.
Expiration Date	Please enquire.



Human ACHE
PicoKine
ELISA Kit
Standard C...