


Product Datasheet

Human C3a anaphylatoxin ELISA Kit (orb865345)

Description	Human C3a anaphylatoxin ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S672-R748
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA, citrate) and urine
Uniprot ID	P01024
Sensitivity	20 pg/ml
Dilution Range	Dilution ratio of 1:2000000, concentration in serum and plasma is 936-1300 ug/ml.
Expiration Date	Please enquire.



Human C3a anaphylatoxin
PicoKine
ELISA K...