

Product Datasheet

heavy chain Myosin/Myh3 Antibody (orb865435)

Description

heavy chain Myosin/Myh3 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, FC, IHC, WB

Applications
Immunogen

E.coli-derived rat heavy chain Myosin/Myh3 recombinant protein (Position: K44-S119).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW


224 kDa

Uniprot ID
[P12847](#)
Dilution Range


Western blot, 0.25-0.5 µg/ml, Rat Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Mouse, Rat Flow Cytometry, 1-3 µg/1x10⁶ cells, Human, Rat Direct ELISA, 0.1-0.5 µg/ml, Rat

Expiration Date


12 months from date of receipt.



Western blot analysis of heavy chain Myo...



IHC analysis of heavy chain Myosin/Myh3 ...



IHC analysis of heavy chain Myosin/Myh3 ...