

Product Datasheet

Rpl9 Antibody (orb865565)

Description

Rpl9 Antibody

Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, FC, ICC, IF, IHC, WB

Immunogen

E.coli-derived mouse Rpl9 recombinant protein (Position: M1-E192).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

22 kDa

Uniprot ID


[P51410](#)

Dilution Range


Western blot, 0.1-0.25 µg/ml, Mouse, Rat
Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml,
Mouse Immunocytochemistry/Immunofluorescence, 5 µg/ml,
Mouse Immunofluorescence, 5 µg/ml, Mouse Flow Cytometry,
1-3 µg/1x10⁶ cells, Mouse Direct ELISA, 0.1-0.5 µg/ml, Mouse

Expiration Date


12 months from date of receipt.



Western blot analysis of Rpl9 using anti-Rpl9 antibody.



IHC analysis of Rpl9 using anti-Rpl9 antibody.



IF analysis of Rpl9 using anti-Rpl9 antibody.