

Product Datasheet

DDX4/MVH Antibody (monoclonal, 3C3) (orb865594)

Description

DDX4/MVH Antibody (monoclonal, 3C3)

Species/Host

Mouse

Reactivity

Mouse, Rat

Conjugation

Unconjugated

Tested Applications

IF, IHC, WB

Immunogen

E.coli-derived human DDX4/MVH recombinant protein (Position: D3-D666).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Mouse IgG2b

Clonality

Monoclonal

MW

79 kDa

Uniprot ID


[Q9NQ10](#)

Dilution Range


Western blot, 0.25-0.5 µg/ml, Mouse, Rat Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Mouse, Rat Immunofluorescence, 5 µg/ml, Mouse, Rat

Expiration Date


12 months from date of receipt.



Western blot analysis of DDX4/MVH using ...



IHC analysis of DDX4/MVH using anti-DDX4...



IHC analysis of DDX4/MVH using anti-DDX4...