

Product Datasheet

**splicing factor 1 Antibody (monoclonal, 7D9E3)
(orb865656)**

Description

splicing factor 1 Antibody (monoclonal, 7D9E3)

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, ICC, IF, IHC, WB

Applications
Immunogen

E. coli-derived human splicing factor 1 recombinant protein (Position: R160-Q266).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Mouse IgG2b

Clonality

Monoclonal

MW


68 kDa

Uniprot ID
[Q15637](https://www.uniprot.org/entry/Q15637)
Dilution Range


 Western blot, 0.25-0.5 µg/ml, Human, Mouse, Rat
 Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human, Mouse, Rat
 Immunocytochemistry/Immunofluorescence, 5 µg/ml, Human Flow Cytometry, 1-3 µg/1x10⁶ cells, Human, Mouse, Rat

Expiration Date


12 months from date of receipt.



Western blot analysis of splicing factor...



IHC analysis of splicing factor 1 using ...



IHC analysis of splicing factor 1 using ...