

Product Datasheet

Filamin B/FLNB Antibody (monoclonal, 11E2D2) (orb865673)

Description

Filamin B/FLNB Antibody (monoclonal, 11E2D2)

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

ICC, IF, IHC, WB

Immunogen

E.coli-derived human Filamin B/FLNB recombinant protein (Position: Q397-D701).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Mouse IgG1

Clonality

Monoclonal

MW

278 kDa

Uniprot ID


[O75369](#)

Dilution Range


Western blot, 0.25-0.5 µg/ml, Human
Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human
Immunocytochemistry/Immunofluorescence, 5 µg/ml, Human

Expiration Date


12 months from date of receipt.



Western blot analysis of Filamin B/FLNB ...



IHC analysis of Filamin B/FLNB using ant...



IHC analysis of Filamin B/FLNB using ant...