


## Product Datasheet

### Hemoglobin/HBA1/HBA2 Antibody (orb865699)


## Description

Hemoglobin/HBA1/HBA2 Antibody


|                            |                                                                                                                                                                                                                                        |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                                                 |
| <b>Reactivity</b>          | Human, Mouse                                                                                                                                                                                                                           |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                           |
| <b>Tested Applications</b> | IHC, WB                                                                                                                                                                                                                                |
| <b>Immunogen</b>           | A synthetic peptide corresponding to a sequence of human Hemoglobin/HBA1/HBA2 (AVHASLDKFLASVSTVLTSKYR).                                                                                                                                |
| <b>Form/Appearance</b>     | Lyophilized                                                                                                                                                                                                                            |
| <b>Concentration</b>       | Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.                                                                                                                                                              |
| <b>Storage</b>             | At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.                                      |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                  |
| <b>Application notes</b>   | Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml. |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                                                                                             |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                             |
| <b>MW</b>                  | 15 kDa                                                                                                                                                                                                                                 |
| <b>Uniprot ID</b>          | <b>P69905</b>                                                                                                                                                                                                                          |
| <b>Dilution Range</b>      | Western blot, 0.1-0.25 µg/ml, Human, Mouse<br>Immunohistochemistry(Paraffin-embedded Section), 1-2 µg/ml, Human                                                                                                                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                        |



Western blot analysis of Hemoglobin/HBA1...



IHC analysis of Hemoglobin/HBA1/HBA2 usi...



IHC analysis of Hemoglobin/HBA1/HBA2 usi...