


Product Datasheet

Myc-Tag Mouse Monoclonal Antibody (orb865746)

Description	Myc-Tag Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IP, WB
Immunogen	A synthesized peptide
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB 1:1000-1:5000 ICC/IF 1:100-1:200 IP 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of Myc-tagged prot...