



Product Datasheet

beta Actin Mouse Monoclonal Antibody(HRP Conjugated) (orb865756)

Description	beta Actin Mouse Monoclonal Antibody(HRP Conjugated)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthesized peptide derived from human beta Actin
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P60709
Dilution Range	WB 1:3000-1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of beta Actin expr...



Western blot analysis of beta Actin expr...