



Product Datasheet

p57 Kip2 Rabbit Monoclonal Antibody (orb865763)


Description	p57 Kip2 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	A synthesized peptide derived from human p57 Kip2
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P49918
Dilution Range	WB: 1:500-1:2000, IHC: 1:50-1:200, ICC/IF 1:50-1:200, IP 1:50
Expiration Date	12 months from date of receipt.



Immunofluorescent analysis of HeLa cells...



Immunohistochemical analysis of paraffin...



Western blot analysis on HeLa cell using...