


Product Datasheet

APOC3 Rabbit Monoclonal Antibody (orb865791)

Description	APOC3 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthesized peptide derived from human APOC3
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P02656
Dilution Range	WB 1:500-1:2000IHC 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of APOC3 expression...