


Product Datasheet

OGG1 Rabbit Monoclonal Antibody (orb865796)

Description	OGG1 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthesized peptide derived from human OGG1
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	O15527
Dilution Range	WB 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of OGG1 expression...