


## Product Datasheet

# Transferrin Receptor 2 Rabbit Monoclonal Antibody (orb865951)

|                            |                                                                                                                                            |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Transferrin Receptor 2 Rabbit Monoclonal Antibody                                                                                          |
| <b>Species/Host</b>        | Rabbit                                                                                                                                     |
| <b>Reactivity</b>          | Human                                                                                                                                      |
| <b>Conjugation</b>         | Unconjugated                                                                                                                               |
| <b>Tested Applications</b> | IHC, WB                                                                                                                                    |
| <b>Immunogen</b>           | A synthesized peptide derived from human Transferrin Receptor 2                                                                            |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                     |
| <b>Concentration</b>       | Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.                                               |
| <b>Storage</b>             | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                      |
| <b>Isotype</b>             | IgG                                                                                                                                        |
| <b>Clonality</b>           | Monoclonal                                                                                                                                 |
| <b>MW</b>                  | 100 kDa                                                                                                                                    |
| <b>Uniprot ID</b>          | <b>Q9UP52</b>                                                                                                                              |
| <b>Dilution Range</b>      | WB 1:500-1:2000IHC 1:50-1:200                                                                                                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                            |



Western blot analysis of Transferrin Rec...