


Product Datasheet

MIB1 / DIP1 Rabbit Monoclonal Antibody (orb865982)

Description	MIB1 / DIP1 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthesized peptide derived from human MIB1 / DIP1
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	110 kDa
Uniprot ID	Q86YT6
Dilution Range	WB 1:500-1:2000IHC 1:50-1:200ICC/IF 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of MIB1/DIP1 expression...