


Product Datasheet

RENT1 / hUPF1 Rabbit Monoclonal Antibody (orb866024)

Description	RENT1 / hUPF1 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, IP, WB
Immunogen	A synthesized peptide derived from human RENT1 / hUPF1
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	130 kDa
Uniprot ID	Q92900
Dilution Range	WB 1:1000-1:5000IHC 1:50-1:200ICC/IF 1:50-1:200IP 1:50FC 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of RENT1/hUPF1 exp...