


Product Datasheet

TNNI2 Rabbit Monoclonal Antibody (orb866157)

Description	TNNI2 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthesized peptide derived from human TNNI2
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	21 kDa
Uniprot ID	P48788
Dilution Range	WB 1:1000-1:5000IP 1:30
Expiration Date	12 months from date of receipt.



Western blot analysis of Troponin I fast...