


Product Datasheet

Nogo Receptor Rabbit Monoclonal Antibody (orb866309)

Description	Nogo Receptor Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthesized peptide derived from human Nogo Receptor
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	66 kDa
Uniprot ID	Q9BZR6
Dilution Range	WB 1:1000-1:5000IP 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of Nogo Receptor e...