


Product Datasheet

Histone H3 (acetyl K18) Rabbit Monoclonal Antibody (orb866688)

Description	Histone H3 (acetyl K18) Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	A synthesized peptide derived from Histone H3 (acetyl K18)
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	17 kDa
Uniprot ID	P68431
Dilution Range	WB 1:500-1:2000IHC 1:50-1:200ICC/IF 1:50-1:200FC 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of Histone H3 (ace...