


Product Datasheet

Histone H3 (acetyl K23) Rabbit Monoclonal Antibody (orb866689)

Description	Histone H3 (acetyl K23) Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ChIP, ICC, IF, WB
Immunogen	A synthesized peptide derived from Histone H3 (acetyl K23)
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	15 kDa
Uniprot ID	Q71DI3, Q6NXT2, P68431, Q16695, P84243
Dilution Range	WB 1:500-1:2000 ICC/IF 1:50-1:200 ChIP 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of Histone H3 (acetyl K23)