

Product Datasheet

ACE2 antibody (orb866860)

Description

Goat polyclonal antibody to ACE2

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FC, IHC, WB

Applications
Target

ACE2 (N Terminal)

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

WB: Approx 90-100kDa band observed in lysates of cell line HepG2 (calculated MW of 92.4kDa according to NP_001358344.1). Recommended concentration 0.5-2ug/ml. Primary incubation 1 hour at room temperature. Approx 100kDa band was also observed on full-lengt. IHC: Paraffin embedded Human Testis. Recommended concentration: 8µg/ml.. Notes: Flow Cytometry: Flow cytometric analysis of HepG2 cells.. Recommended concentration: 10ug/ml.

Protein Sequence

PROPRIETARY,

Clonality

Polyclonal

MW

92.4

Entrez
59272
Dilution Range

Peptide ELISA: antibody detection limit dilution 1:128000.

Expiration Date

12 months from date of receipt.