

Product Datasheet

ACE2 antibody (orb866862)

Description	Goat polyclonal antibody to ACE2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	ACE2 (N Terminal)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: Paraffin embedded Human Testis. Recommended concentration: 8µg/ml. A customer reported weak staining on paraffin embedded Human Kidney, at a concentration of 5ug/ml.
Protein Sequence	PROPRIETARY
Clonality	Polyclonal
MW	92.4
Entrez	59272
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000.
Expiration Date	12 months from date of receipt.