

## Product Datasheet

### Human DES ELISA Kit (orb866889)

|                            |                                                              |
|----------------------------|--------------------------------------------------------------|
| <b>Description</b>         | Sandwich Human DES ELISA Kit                                 |
| <b>Reactivity</b>          | Human                                                        |
| <b>Range</b>               | 0.16-10 ng/mL                                                |
| <b>Tested Applications</b> | ELISA                                                        |
| <b>Note</b>                | For research use only                                        |
| <b>Sample Types</b>        | Tissue homogenates, cell lysates and other biological fluids |
| <b>Sensitivity</b>         | 0.062 ng/mL                                                  |
| <b>Expiration Date</b>     | Please enquire.                                              |

