

## Product Datasheet

### Canine CHGA ELISA Kit (orb866929)

<b>Description</b>	Sandwich Canine CHGA ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	31.25-2000 pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, tissue homogenates and other biological fluids
<b>Sensitivity</b>	57 pg/mL
<b>Expiration Date</b>	Please enquire.

