

## Product Datasheet

### Fish GnRH ELISA Kit (orb866958)

|                            |                                                                                                      |
|----------------------------|------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Sandwich Fish GnRH ELISA Kit                                                                         |
| <b>Reactivity</b>          | Fish                                                                                                 |
| <b>Range</b>               | 0.25-16 ng/mL                                                                                        |
| <b>Tested Applications</b> | ELISA                                                                                                |
| <b>Note</b>                | For research use only                                                                                |
| <b>Sample Types</b>        | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| <b>Sensitivity</b>         | 0.1 ng/mL                                                                                            |
| <b>Expiration Date</b>     | Please enquire.                                                                                      |

