

Product Datasheet

SARS-CoV-2 SPIKE PROTEIN S1 antibody (orb867040)

Description	Rabbit polyclonal antibody to SARS-CoV-2 SPIKE PROTEIN S1
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC-P, WB
Immunogen	Peptide corresponding to 16 amino acids near the amino terminus of SARS-CoV-2 (COVID-19) Spike S1 glycoprotein. The immunogen is located within the first 50 amino acids of SARS-CoV-2 (COVID-19) Spike S1 protein.
Target	SARS-CoV-2 SPIKE PROTEIN S1
Preservatives	0.02% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.