

Product Datasheet

MMP-3 antibody (Biotin) (orb867066)

Description	Mouse monoclonal antibody against MMP-3 conjugated to Biotin
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	E. coli derived recombinant human partial length of MMP-3 (Phe100-Cys477)
Target	MMP-3
Preservatives	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	IgG concentration 0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.