

Product Datasheet

Human INTERLEUKIN-11 protein (orb867201)

Description	Recombinant human INTERLEUKIN-11 protein
Reactivity	Human
Tested Applications	ELISA
Target	IL-11
Preservatives	1.0% Trehalose0.05% Tween80. 50 mM Tris-HCl, pH 8.0, 0.15 M NaCl
Form/Appearance	Purified recombinant protein - lyophilized
Concentration	100 µg/ml after reconstitution
Storage	-70°C
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.