

## Product Datasheet

### Human CCL24 protein (orb867203)

|                            |                                                                |
|----------------------------|----------------------------------------------------------------|
| <b>Description</b>         | Recombinant human CCL24 protein                                |
| <b>Reactivity</b>          | Human                                                          |
| <b>Tested Applications</b> | ELISA                                                          |
| <b>Target</b>              | CCL24                                                          |
| <b>Preservatives</b>       | 1.0% Trehalose. 50 mM Tris-HCl, pH 8.0, 0.15 M Sodium Chloride |
| <b>Form/Appearance</b>     | Purified recombinant protein - lyophilized                     |
| <b>Concentration</b>       | 100 µg/ml after reconstitution                                 |
| <b>Storage</b>             | -70°C                                                          |
| <b>Note</b>                | For research use only                                          |
| <b>Clonality</b>           | Recombinant                                                    |
| <b>Purity</b>              | Rec. Protein                                                   |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                 |