

Product Datasheet

Human CCL17 protein (orb867204)

Description	Recombinant human CCL17 protein
Reactivity	Human
Tested Applications	ELISA
Target	CCL17
Preservatives	1.0% Trehalose. 50 mM Tris-HCl, pH 8.0, 0.15 M Sodium Chloride
Form/Appearance	Purified recombinant protein - lyophilized
Concentration	100 µg/ml after reconstitution
Storage	-70°C
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.