

Product Datasheet

PIN-Cocktail 2B-P504S antibody [Out of stock] (orb867255)

Description	Mouse monoclonal antibody to PIN-Cocktail 2B-P504S
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P
Immunogen	1. Recombinant human AMACR (P504S) protein, whole sequence; 2. Recombinant human p63 protein aa 1-203; 3. Cytokeratin extract of human stratum corneum
Preservatives	Affinity purified antibodies in PBS, rAlbumin, sodium azide (0.09%)**. Use antibody dilution buffer (e.g. Art. No. PU002) containing sufficient protein and preservative.
Storage	2-8°C
Note	For research use only
Isotype	IgG/IgG1
Dilution Range	(liquid conc.) 1:25-1:50
Expiration Date	12 months from date of receipt.