


Product Datasheet

Mouse IL-15 R alpha / CD215 Protein (orb867324)

Description

Mouse IL-15 R alpha, Fc Tag is expressed from

	Lys ...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-15 R alpha
Preservatives	50 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 44.8 kDa. The protein migrates as 60-90 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	Mouse IL-15 R alpha, Fc Tag (orb867324) is expressed from human 293 cells (HEK293). It contains AA Gly 33 - Lys 205 (Accession # Q60819-1).
MW	44.8 kDa
Uniprot ID	Q60819
NCBI	NP_032384.1
Expiration Date	6 months from date of receipt.



Immobilized
Human IL-
15,
premium
grade (...)