

Product Datasheet

Human Latent TGF-beta 2 / TGFB2 Protein (orb867336)

Description

Human Latent TGFB2, His Tag (TG2-

cells (HEK293). It contains AA L...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Latent TGF-beta 2

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

N-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 47.6 kDa. The protein migrates as 12 kDa and 40-50 kDa under reducing (R) condition (SDS-PAGE).

Protein Sequence

NP_003229.1

Purity

● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Source

Human Latent TGFB2, His Tag (orb867336) is expressed from human 293 cells (HEK293). It contains AA Leu 21 - Ser 414 (Accession # P61812-1).

MW

47.6 kDa

Uniprot ID
P61812
NCBI
NP_003229.1
Expiration Date

6 months from date of receipt.