

Product Datasheet

Mouse CD39 / ENTPD2 Protein (orb867339)

Description	Mouse CD39, His Tag is expressed from human 293 cells (HEK293). It contains AA Thr 38 -Ile 478 (Acce...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CD39
Preservatives	20mM Tris, 150mM NaCl, pH8.0
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 51.7 kDa. The protein migrates as 66-70 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_033978.1
Purity	95%
Source	Mouse CD39, His Tag (orb867339) is expressed from human 293 cells (HEK293). It contains AA Thr 38 -Ile 478 (Accession # P55772-1).
MW	51.7 kDa
Uniprot ID	P55772
NCBI	NP_033978.1
Expiration Date	6 months from date of receipt.