

Product Datasheet

Human Tau protein (orb867352)

Description	Recombinant Human Tau protein
Endotoxins	Less than 0.1 EU/ug
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-His
Note	For research use only
Protein Sequence	Ala2-Leu441
Source	Human 293 Cells
Uniprot ID	P10636
Expiration Date	12 months from date of receipt.