

Product Datasheet

Human SIGLEC-8 protein (orb867414)

Description	Recombinant Human SIGLEC-8 protein
Endotoxins	Less than 1 EU/ug
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Human Fc
Note	For research use only
Protein Sequence	Met17-Ala363
Source	Human Cells
Uniprot ID	Q9NYZ4
Expiration Date	12 months from date of receipt.