

Product Datasheet

Anti-NKCC1/SLC12A2/PE-Cy5.5 (orb870407)

Description	Anti-NKCC1/SLC12A2/PE-Cy5.5
Species/Host	Rabbit
Reactivity	Human
Conjugation	PE/Cy5.5
Tested Applications	FC, ICC
Immunogen	KLH conjugated synthetic peptide derived from human NKCC1
Target	NKCC1/SLC12A2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P55011
Dilution Range	Flow-Cyt=3ug/test,ICC=1:100-500
Expiration Date	12 months from date of receipt.