

Product Datasheet

Anti-CD73/PE-Cy7 (orb873792)

Description	Anti-CD73/PE-Cy7
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	PE/Cy7
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human CD73/NT5
Target	CD73
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q61503
Expiration Date	12 months from date of receipt.