

Product Datasheet

Anti-Phospho-Smad3 (Ser208)/PE-Cy5 (orb879134)

Description	Anti-Phospho-Smad3 (Ser208)/PE-Cy5
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	PE/Cy5
Tested Applications	FC, ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Smad3 around the phosphorylation site of Ser208
Target	Phospho-Smad3 (Ser208)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P84022
Dilution Range	Flow-Cyt=1µg/Test,ICC=1:100
Expiration Date	12 months from date of receipt.