

Product Datasheet

Myelin Basic Protein antibody [Out of stock] (orb88041)

Description Mouse Monoclonal to Myelin Basic Protein.

Species/Host Mouse

Reactivity Goat, Human, Rabbit, Rat, Sheep

Conjugation Unconjugated

Tested Applications ELISA, IHC-Fr, IHC-P, WB

Immunogen Purified human myelin basic protein.

Target MBP

Preservatives PBS, pH 7.4, no preservatives.

Concentration 0.96 mg/ml

Storage Short term: store at 4°C. Long term: aliquot and store at -20°C. Avoid freeze-thaw cycles.

Note For research use only


Isotype IgG2b

Clonality Monoclonal

Entrez [4155](#)

Dilution Range ELISA:10-20 µg/ml,IHC-P:5 µg/ml,WB:1-2 µg/ml

Expiration Date 12 months from date of receipt.



Immunohistochemical staining of paraffin...