

Product Datasheet

DARPP32 (phospho-Thr34) antibody (FITC) [Out of stock] (orb8808)

Description

Rabbit polyclonal to DARPP32

conjugated to FITC

Species/Host

Rabbit

Conjugation

FITC

Tested

IF

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human DARPP32 around the phosphorylation site of Thr34

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

32kDa

Uniprot ID
[Q9UD71](#), [Q60829](#), [Q6J4I0](#)
Dilution Range

IF: 1:50-200

Expiration Date

12 months from date of receipt.