

Product Datasheet

Anti-human Fibrinogen/PE-Cy5 (orb883201)

Description	Anti-human Fibrinogen/PE-Cy5
Species/Host	Mouse
Reactivity	Human
Conjugation	PE/Cy5
Tested Applications	FC
Immunogen	Fibrinogen from human plasma
Target	human Fibrinogen
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
CAS Number	9001-32-5
Uniprot ID	P02671
Expiration Date	12 months from date of receipt.