

Product Datasheet

Anti-Phospho-PKR (Thr446 + Thr451)/PE-Cy5 (orb886925)

Description

Anti-Phospho-PKR (Thr446 +

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE/Cy5
Tested Applications	FC, ICC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human PKR around the phosphorylation site of Thr446/451
Target	EIF2AK2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	62kDa
Uniprot ID	P19525
Dilution Range	Flow-Cyt=2ug/Test,ICC=1:100
Expiration Date	12 months from date of receipt.