

Product Datasheet

**Anti-phospho-PI3 kinase p85 alpha + gamma
(Tyr467 + Tyr199)/PE-Cy5 (orb887079)**

Description

Anti-phospho-PI3 kinase p85

Species/Host

Tyr199)/PE-Cy5

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

PE/Cy5

Tested

FC

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human PI3 kinase p85 alpha around the phosphorylation site of Tyr467

Target

phospho-PI3 kinase p85 alpha + gamma (Tyr467 + Tyr199)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID

P27986

Expiration Date

12 months from date of receipt.