

Product Datasheet

Anti-cGKI/AP (orb887315)

Description

Anti-cGKI/AP

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	AP
Tested Applications	ICC
Immunogen	KLH conjugated synthetic peptide derived from human Prkg1/cGKI corresponding to amino acids 582-628.
Target	PRKG1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	74kDa
Uniprot ID	Q13976
Dilution Range	WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Expiration Date	12 months from date of receipt.