

## Product Datasheet

### phospho-c-Jun (Ser243)/PE-Cy5 antibody [Out of stock] (orb888497)

|                            |                                                                                                              |
|----------------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | phospho-c-Jun (Ser243)/PE-Cy5 antibody                                                                       |
| <b>Species/Host</b>        | Rabbit                                                                                                       |
| <b>Conjugation</b>         | PE/Cy5                                                                                                       |
| <b>Tested Applications</b> | FC                                                                                                           |
| <b>Immunogen</b>           | KLH conjugated Synthesised phosphopeptide derived from human c-Jun around the phosphorylation site of Ser243 |
| <b>Target</b>              | phospho-c-Jun (Ser243)                                                                                       |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.                                        |
| <b>Form/Appearance</b>     | Liquid                                                                                                       |
| <b>Concentration</b>       | 1mg/ml                                                                                                       |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.                             |
| <b>Note</b>                | For research use only                                                                                        |
| <b>Isotype</b>             | IgG                                                                                                          |
| <b>Clonality</b>           | Polyclonal                                                                                                   |
| <b>MW</b>                  | 43kDa                                                                                                        |
| <b>Uniprot ID</b>          | <b>P05412</b>                                                                                                |
| <b>Dilution Range</b>      | Flow-Cyt=0.2µg /Test                                                                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                              |