

Product Datasheet

Phospho-FAK (Tyr397)/PE-Cy5 antibody [Out of stock] (orb888674)

Description	Phospho-FAK (Tyr397)/PE-Cy5 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE/Cy5
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human FAK around the phosphorylation site of Tyr397
Target	Phospho-FAK (Tyr397)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	116kDa
Uniprot ID	Q05397
Expiration Date	12 months from date of receipt.