

Product Datasheet

PLC gamma 1 (phospho-Tyr771) antibody (FITC)
[Out of stock] (orb8889)

Description

Rabbit polyclonal to PLC gamma 1

Species/Host

Rabbit

Conjugation

FITC

Tested

IF

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human PLC gamma 1 around the phosphorylation site of Tyr771

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

148kDa

Uniprot ID
P10686, Q62077, P19174
Dilution Range

IF: 1:50-200

Expiration Date

12 months from date of receipt.