

Product Datasheet

Anti-Beta tubulin /AP (orb891753)

| | |
|----------------------------|--|
| Description | Anti-Beta tubulin /AP |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | AP |
| Tested Applications | IHC-P, WB |
| Immunogen | KLH conjugated synthetic peptide derived from human tubulin Beta |
| Target | Beta tubulin (Loading Control) |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P07437 |
| Dilution Range | WB=1:5000-20000,IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000 |
| Expiration Date | 12 months from date of receipt. |